

WP 300.11

Device for cupping tests



Description

 investigation of the material properties of thin sheet metal
cupping test based on DIN 50101 T1

This accessory for WP 300 makes it possible to investigate the material properties of thin sheet metal. To do this, the specimens are deformed in practice-oriented cupping experiments and then examined.

The accessory contains a punch, a holddown clamp, a die and a mirrored cylinder for analysing the specimen. The assembled test device is installed in the compression area of the WP 300 experimental unit, between bottom crossbar and crosshead.

A thin sheet metal specimen is clamped between the die and the hold-down device. The specimen in then deformed until it cracks by means of a spherical drawing punch. The mirror allows you to observe the specimen from below and to see exactly when the crack starts. Thin sheet metal with a thickness up to 2mm can be studied with the device.

The depth achieved in mm is defined as a standard of comparison for cold formability. In addition, the type of the rupture and the surface structure of the sagging area are analysed.

A set of specimens is included to carry out experiments. The set contains five specimens each of aluminium, copper, brass and steel.

Learning objectives/experiments

 determine the quality of thin sheet metal made of various materials

Specification

- [1] cupping test based on DIN 50101 T1 for quality assessment of thin sheet metal
- [2] 20 specimens of steel, aluminium, copper, brass
- [3] spherical drawing punch
- [4] mirror for analysing the underside of the specimen
- [5] accessory for WP 300

Technical data

5 specimens each of aluminium, copper, brass, steel

specimen dimensions
60x60x1mm

Drawing punch: sphere 20mm Die diameter: Ø 27mm

LxWxH: 110x60x250mm

Scope of delivery

Weight: approx. 3kg

- 1 device for cupping tests
- 1 set of specimens (20 pieces)
- 1 set of accessories

G.U.N.T. Gerätebau GmbH, Hanskampring 15-17, D-22885 Barsbüttel, Telefon (040) 67 08 54-0, Fax (040) 67 08 54-42, Email sales@gunt.de, Web www.gunt.de We reserve the right to modify our products without any notifications. Page 1/2 - 11.2023



WP 300.11 Device for cupping tests

Required accessories

WP 300 Materials testing, 20kN

Optional accessories

Set of 5 cupping specimens, Al
Set of 5 cupping specimens, Cu
Set of 5 cupping specimens, St
Set of 5 cupping specimens, CuZn

G.U.N.T. Gerätebau GmbH, Hanskampring 15-17, D-22885 Barsbüttel, Telefon (040) 67 08 54-0, Fax (040) 67 08 54-42, Email sales@gunt.de, Web www.gunt.de We reserve the right to modify our products without any notifications. Page 2/2 - 11.2023