


# WP 300.25

Set of 4 tension specimens, flat, Al, Cu, St, CuZn



#### Description

- tensile specimens of rectangular cross-section in accordance with DIN 50125, shape E
- together with WP 300.14 as accessories for the WP 300 experimental unit

This set of specimens is available as an accessory for the WP 310 experimental unit. The specimens comply with the requirements of DIN 50125, shape E and are of rectangular cross-section. For marking the test length, there are two centre punch marks on the shafts.

The specimen is installed in the clamping device of WP 300.14. The specimens have a widened section at both ends for attachment, using this section the specimens can be hung in the clamping device. The resulting device is clamped into the tensile area of WP 310, between top crossbar and crosshead. In the tensile test, a uniaxial state of stress is produced in the specimen by an external load in the longitudinal direction. Then a uniform normal stress distribution prevails in the test cross-section of the specimen. In order to determine the strength of the material, the load on the specimen is slowly and steadily increased, until the specimen breaks.

The set contains one specimen made of aluminium, one of copper, one of brass and one of steel.

### Learning objectives/experiments

 together with WP 300.14: tensile tests with flat specimens of different materials on the WP 300 experimental unit

## Specification

- [1] flat specimens for tensile tests
- [2] 4 specimens madee of aluminium, copper, brass, steel
- [3] specimens according to DIN 50125, shape E
- [4] accessory for clamping device WP 300.14
- [5] accessory for WP 300 Materials testing, 20kN

#### **Technical data**

4 tensile specimens in acc. with DIN 50125, shape E

- LxWxH: 144x45x1,5mm
- cross-section: 12,5x1,5mm
- test length: 25mm
- total specimen length: 144mm
- material: Al, Cu, St, CuZn

Weight: approx. 100g

### Scope of delivery

1 set of specimens (4 pieces)

G.U.N.T. Gerätebau GmbH, Hanskampring 15-17, D-22885 Barsbüttel, Telefon (040) 67 08 54-0, Fax (040) 67 08 54-42, Email sales@gunt.de, Web www.gunt.de We reserve the right to modify our products without any notifications. Page 1/2 - 11.2023



# WP 300.25 Set of 4 tension specimens, flat, Al, Cu, St, CuZn

Required accessories

WP 300Materials testing, 20kNWP 300.14Clamping device for flat tensile specimens

G.U.N.T. Gerätebau GmbH, Hanskampring 15-17, D-22885 Barsbüttel, Telefon (040) 67 08 54-0, Fax (040) 67 08 54-42, Email sales@gunt.de, Web www.gunt.de We reserve the right to modify our products without any notifications. Page 2/2 - 11.2023